

Nira

Garanti: 12 ay
Bu ¼r¼ne ait T¼rkiye i¼i kapı teslim fiyat ve teslim s¼resi i¼eren teklifimizi almak i¼in **info@yursat.com.tr** e-posta adresine baŐvuru yapabilir ya da ayrıntılı bilgi i¼in +90 224 240 03 04 numaralı telefonumuzdan bizlere ulaŐabilirsiniz.
Nira Markası, tedarik s¼resi i¼in l¼tfen bizimle iletiŐime ge¼iniz.

Firmamız Nira T¼rkiye Distrib¼t¼r¼ veya temsilcisi deĐildir. Firmamız sipariŐ durumunda, belirtilen ¼r¼nlerde sadece Orjinal ve yeni ¼r¼n teklifi sunmaktadır. Bu sitede g¼sterilen ¼zel marka adları ve ticari markalar ilgili sahiplerinin m¼lkiyetindedir, talep durumunda kaldırılmaktadır.



¼r¼n	A¼ıklama
Neptun 800	
Neptune 803 Table Top Version	Neptune 803 comes in two layouts, table top and process version. The core of the instruments and their analytical characteristics are exactly the same. Table top version, the most compact version of Neptune803, is designed to be installed in laboratories or production areas and is complete with: - Capillary column 50m x0,53mm. - Connection kit, Air and hydrogen pipes(5m) and cables(2m). - Laser Printer, USB. - Set of n°50 vials, Pyrex® glass. - Syringe for standard mixture injection, 5µl. - Standard mixture (up to 15 solvents) for calibration 5ml. - First calibration at our laboratories with your std. mixture of your choice. - Package: Wood Box. Note: this equipment shall be installed in an environmental temperature range between +10 and +35°C. Technical characteristics: Detector : flame ionization Desorber (Headspace) : accepts 20 ml sample receptacles Incubation times : 0, 5, 10, 30 minutes Working temperature : +10, +40°C(*) Desorber temperature : from 50 to 150°C
Analyzer Mod. Neptune 801 Table-Top Version	

[Analyzer Mod. Neptune 801 Process
Version](#)