

Nira

Garanti: 12 ay
Bu rne ait Trkiye ii kapı teslim fiyat ve teslim sresi ieren teklifimizi almak iin **info@yursat.com.tr** e-posta adresine baŐvuru yapabilir ya da ayrıntılı bilgi iin +90 224 240 03 04 numaralı telefonumuzdan bizlere ulaŐabilirsiniz.
Nira Markası, tedarik sresi iin ltfen bizimle iletiŐime geiniz.

Firmamız Nira Trkiye Distribtr veya temsilcisi deĐildir. Firmamız sipariŐ durumunda, belirtilen rnlerde sadece Orjinal ve yeni rn teklifi sunmaktadır. Bu sitede gsterilen zel marka adları ve ticari markalar ilgili sahiplerinin mlkiyetindedir, talep durumunda kaldırılmaktadır.



| rn | Aıklama |
|---|---|
| Neptun 800 | |
| Neptune 803 Table Top Version | Neptune 803 comes in two layouts, table top and process version. The core of the instruments and their analytical characteristics are exactly the same. Table top version, the most compact version of Neptune803, is designed to be installed in laboratories or production areas and is complete with: - Capillary column 50m x0,53mm. - Connection kit, Air and hydrogen pipes(5m) and cables(2m). - Laser Printer, USB. - Set of n50 vials, Pyrex glass. - Syringe for standard mixture injection, 5l. - Standard mixture (up to 15 solvents) for calibration 5ml. - First calibration at our laboratories with your std. mixture of your choice. - Package: Wood Box. Note: this equipment shall be installed in an environmental temperature range between +10 and +35C. Technical characteristics: Detector : flame ionization Desorber (Headspace) : accepts 20 ml sample receptacles Incubation times : 0, 5, 10, 30 minutes Working temperature : +10, +40C(*) Desorber temperature : from 50 to 150C |
| Analyzer Mod. Neptune 801 Table-Top Version | |

[Analyzer Mod. Neptune 801 Process
Version](#)